USER MANUAL

Lumbosacral corset with supporting straps
Thoracolumbosacral corset with supporting straps
Product code: LumboTim/LumboTim Plus

1. Characteristics

The corset is made of elastic compression fabric fastened by velcro at the front. It consists of flexible underwires and aluminium stays, which relieve and stabilize the spine, and a tightening system which adjusts the force of stabilization. Owing to the spatial structure of compression fabric, the corset is air permeable and removes excessive heat and moisture perfectly. The height of the corset is appropriate to stabilize the lumbar spine (LumboTim) or the thoracic-lumbosacral spine (LumboTim Plus).

2. Application

The corset is used to reduce pain, in decompression of nerve roots, in local stabilization of spine, as well as in strengthening during verticalization of a patient who suffers from orthopedic or neurologic injury. Owing to elastic fabric, the material has compression effect on the muscles within the trunk.

Indications for use:

- Discopathy
- Spondylolisthesis
- Inflammations of spine joints
- Overload pain
- Degenerative changes
- Neuropathic and radiculopathy
- Increased and painful states of back muscle tone
- In post-surgery treatment and rehabilitation process

3. Contraindications

Do not use in case of:

- Skin abrasion
- Eczema
- Wound and other skin injuries
- After application of warming gel or rub
- Individual contraindications which were determined by the doctor during interview and physical examination

4. Notes

Timago International Group does not bear responsibility of improper use of the product, failure to comply safety regulations and misuse.

5. Fitting instruction

Note! Follow the information contained in this manual while fitting and adjustment. Users with severe pain in the lumbar spine and with reduced mobility may need help from a third person who has appropriate skills.

- To put the corset on, unfasten the velcro and supporting straps which are located at the front of the corset. Wrap the trunk around from the back to the front at the level of the lumbosacral spine. Fasten the stabilizing strap at the front and lateral supporting straps.
- Remember that only properly matched size guarantees correct stabilization and relief.

6. Technical data

SIZE	TRUNK CIRCUMFERENCE (cm)
S	75-85
M	85-95

L	95-105
XL	105-115
XXL	115-125
XXXL	125-140

7. Maintenance

- Hand wash in warm water at 30°C. Do not: tumble dry, chlorinate, iron, bleach, dry clean.
- Ensure the velcro is fastened in order not to destroy the material.
- Dry at room temperature far from intensive heat source.

Note! Maintenance information can be found on the material label.

8. Storage and transport conditions

Product should not be exposed to water or moisture.

9. Warranty information

All products distributed by our company are covered by a warranty, the terms of which are described in the warranty card available on our website. Please contact the retailer from whom you purchased the product. Please, keep in mind that for warranty purposes, the proof of purchase (receipt or invoice) should be retained.

10. Labels

REF	Reference number		Manufacturer
LOT	LOT number	UDI	UDI code
SN	Serial number	MD	Medical device

	Note		Please read the instruction
CE	The manufacturer checked for compliance with essential medical devices requirements.		

Our complete line of products is available from our distributors. Learn more at: **www.timago.com**
Thank you for choosing Timago!

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