INSTRUCTIONS FOR USE

Elastic wrist support

Product code: TGO-C OSN 407

1. Characteristics

Elastic wrist support is made of seamless orthopaedic compression fabric. The thumb support prevents the band from sliding down. Anatomically profiled aluminium stay and elastic strap stabilize the wrist joint from palmar side. Modern and spatial fabric is air permeable and removes excessive heat and moisture limiting the development of bacterial flora. The wrist support can be applied either on the left or the right hand.

2. Application

Wrist support is used to stabilize the joint and soothe pain which occurs in physical strains i.a. in sport.

Indications for use:

- Pain syndrome of joints caused by inflammations, overloads and micro-injuries
- Advanced degenerative changes of radio-carpal joint and wrist joint
- Minor twists and strains
- Chronic instability

Contraindications:

Do not use in case of:

- Skin abrasion
- Eczema
- Wound and other skin injuries
- After application of warming gel or ointment
- Individual contraindications which were determined by the doctor during interview and physical examination

3. Notes

Timago International Group does not bear responsibility of improper maintenance of the product, failure to comply safety regulations and misuse.

4. Fitting instruction

- Unfasten circumferential strap.
- Depending on application side, insert aluminium stay into the pocket located inside the elastic support.
 Example: when the elastic support is used on the left hand, the bending in aluminium stay should be on the inside of the hand.
- Pull the elastic support on the wrist joint and place the thumb into the thumb hole. Stabilize the wrist with use of circumferential strap fastened by Velcro.

Note! Only properly matched size guarantees reduction of pain.
Note! Long-term use may result in reduction of muscle tone and dysfunction of the limb.

5. Technical data

SIZE	CIRCUMFERENCE WRIST (CM)
S	15-17
M	17-19
L	19-21
XL	21-23

6. Maintenance

- Hand wash in warm water. Do not tumble, chlorinate, iron, bleach, dry clean, use tumble dryer.
- Before washing, take out the aluminium stay.
- Dry at room temperature far from intensive heat source.

Note! Maintenance information can be found on the material label.

7. Storage and transport conditions

The product should not be exposed to direct sunlight.

8. Warranty information

All products distributed by our company are covered by a warranty, the terms of which are described in the warranty card available on our website. Please contact the retailer from whom you purchased the product. Please, keep in mind that, for warranty purposes, the proof of purchase (receipt or invoice) should be retained.

9. Labels

	Reference number		Manufacturer
	LOT number		UDI code
	Serial number		Medical device
	Note		Please read the instruction
	The manufacturer checked for compliance with essential medical devices requirements.		

Our complete line of products is available from our distributors. Learn more at: **www.timago.com**. Thank you for choosing Timago!

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