USER MANUAL
Round shower chair – tool free
Product code: TGR-R KP-O 3432

1. Characteristics
The shower chair features an aluminium frame and a round seat made of durable PVC. Thanks to its height-adjustable legs, the chair can be adapted to the user's height. Additionally, it features a Tool-Free (TF) system that facilitates easy assembly.

The set includes: seat, base tubes, height adjustable legs, mounting screws.

Elements of the shower chair:
- **A.** seat
- **B.** base tubes
- **C.** height adjustment – TF system
- **D.** anti-sliding tips

2. Application
The shower chair is ancillary equipment designed to assist users in performing personal hygiene activities. It is recommended for individuals with reduced mobility that prevents or limits independent performance of daily living tasks. The chair can be used in showers, bathtubs, bathrooms, and other humid environments at home or in healthcare facilities.

3. Assembly
Step 1: Place the seat upside down on a flat surface.
Step 2: Attach the crossed base tubes to the seat using the mounting screws and washers.
Step 3: Insert the legs into the base tubes.

HEIGHT ADJUSTMENT
- To adjust the height of the chair press the locking pin, select appropriate height of the leg and then lock it again, repeating steps on the other legs. The legs are adjustable in 2.5 cm increments.

Note! All legs must be set on the same height level.

4. Technical data

Height	36-54 cm
Seat diameter	32 cm
Weight	1,6 kg
Maximum capacity	136kg

5. Notes

Timago International Group does not bear responsibility of improper use of the product, failure to comply safety regulations and misuse.

6. Maintenance

- Clean the shower chair regularly using only water, sponge and mild detergents for plastic materials such as soap. Dry thoroughly.
- **Do not use** for cleaning: strong detergents, abrasive agents, bleach, petrol or solvents (it may cause frame damage, visible defects and in consequence - loss of warranty rights).
- At least once every three months check condition of the tubes, locking screws and tightening of screws.
- It is recommended to check rubber tips regularly. In case of worn-out, replace the tips with the new ones.

7. Storage and transport conditions

The product should not be exposed to direct sunlight, water or moisture.

8. Warranty information

All products distributed by our company are covered by a warranty, the terms of which are described in the warranty card available on our website. Please contact the retailer from whom you purchased the product. Please, keep in mind that for warranty purposes, the proof of purchase (receipt or invoice) should be retained.

9. Labels

	Reference number		Manufacturer

LOT	LOT number	UDI	UDI code
SN	Serial number	MD	Medical device
	Note		Please read the instruction
CE	The manufacturer checked for compliance with essential medical devices requirements.		

Our complete line of products is available from our distributors. Learn more at: **www.timago.com**. Thank you for choosing Timago!

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