USER MANUAL
Aluminium 4-wheeled rollator
Product code: TGR-R RA 885

1. Characteristics
The 4-wheeled rollator is made of an aluminium frame coated in metallic paint. It is equipped with parking brakes, height-adjustable seat and handles. The product comes with shopping bag, backrest and solid wheels with a narrow 8'' cross-section that can be dismantled for transport.

The set includes: frame with seat, backrest, shopping bag, push handles with rear wheels, front wheels**,** locking knobs (2 pcs)**,** locking clip (2 pcs).

Elements of the rollator:

- **A.** Push handle
- **B.** Brake
- **C.** Locking knob
- **D.** Locking pin
- **E.** Backrest
- **F.** Seat
- **G.** Shopping bag
- **H.** Locking pin – height adjustment
- **I.** Locking clip – height adjustment
- **J.** Cane holder
- **K.** Upper adjustment of brake
- **L.** Lower adjustment of brake

2. Application
The rollator is designed for patients with reduced mobility who need stable support and assistance while walking. It is also recommended for the elderly as an aid in maintaining balance, improving walking comfort, and ensuring protection against falls. The use of the rollator has a positive effect on restoring activity in people with mobility impairments and difficulties standing up. It provides improved stabilization, support, and safety while walking.

3. Assembly

UNFOLDING AND FOLDING THE ROLLATOR

- To unfold the rollator, pull the frame apart and press down on the seat until the rollator is fully opened.
- The rear wheels and push handles are connected by the brake cables. First, insert the push handles into the mounting openings at the top of the rollator and secure them with the locking knobs. Then insert the wheels into the openings at the bottom and secure them with the locking pins.
 Note! Make sure all the wheels are on the same height.
- Insert the front wheels into the front part of the frame and secure them with the locking clips.
- Place the backrest on the upper part of the frame and secure it with the locking pins.
- Remove the protective plastic from the wheels.
- To fold the rollator, lift the seat and the bar located underneath it.

BRAKES

- To stop the rollator, pull up the brake handle and hold.
- While stop or using the seat, lock the parking brake by pushing down the brake handle.
- To release the brake, pull the brake handle upwards.

HEIGHT ADJUSTMENT OF HANDLES

- To adjust the height of the handles, unscrew the locking knob located on the outer side of the frame, select appropriate height and tighten the locking knob afterwards.

Note! The handles must be set on the same level.

HEIGHT ADJUSTMENT OF SEAT

- To adjust the height of the seat, remove the locking clip on front frame and press the locking pin on the back frame. Select appropriate height and block the wheels with locking clip and pin.

ADJUSTMENT OF BRAKES

- Unscrew the nut of the cable, pull it up to the appropriate length and lock it by tightening the locking nut.

- Loosen the mounting screw with the use of open-end wrench, pull the black brake cable to adjust the tension of the cable, and then, tighten the mounting screw.

4. Technical data

Overall width	65 cm
Overall depth	69 cm
Seat dimensions	40 x 25,5 cm
Seat height from the ground	54-62 cm
Handles height adjustment	90-100 cm
Weight	6,8 kg
Maximum capacity	130 kg
Maximum capacity of the shopping bag	5 kg
Frame colour	graphite

5. Notes

Timago International Group does not bear responsibility of improper use of the product, failure to comply safety regulations and misuse.

6. Cleaning and maintenance

- Clean the frame with damp cloth and dry thoroughly.
- Do not use strong polishing agents due to the damage of protecting coatings.
- Clean the dirt from moving parts of the rollator regularly (hinges, clamping screws, locks).

7. Storage and transport conditions

The product should not be exposed to direct sunlight, water or moisture

8. Warranty information

All products distributed by our company are covered by a warranty, the terms of which are described in the warranty card available on our website. Please contact the retailer from whom you purchased the product. Please, keep in mind that for warranty purposes, the proof of purchase (receipt or invoice) should be retained.

9. Labels

REF	Reference number		Manufacturer
LOT	LOT number	UDI	UDI code
SN	Serial number	MD	Medical device
	Note		Please read the instruction
CE	The manufacturer checked for compliance with essential medical devices requirements.		

Our complete line of products is available from our distributors. Learn more at: **www.timago.com**. Thank you for choosing Timago!

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